# HANDHELD PULSE OXIMETER

## • Instructions for use


- ES Pulsooxímetro de mano • Instrucciones de uso
- FR Oxymètre de pouls de poche • Mode d'emploi
- DE HS-Einhandpulsoximeter • Gebrauchsanweisung
- IT Pulsossimetro a mano • Istruzioni per l'uso
- PT Pulsoxímetro portátil • Instruções de utilização

# 3x1.5V AA

IPX1

(manufacturer)

**HENRY SCHEIN INC.**  
135 DURYEA ROAD  
Melville, NY 11747 USA

EC REP

**HENRY SCHEIN SERVICES GmbH**  
MONZASTRASSE 2A  
D-63225 LANGEN, GERMANY  
Made in China

Rev.2019/07

## Copyright

Our company owns all rights of this unpublished work and intends to maintain it as a confidential work. We may also seek to maintain this work as an unpublished copyright. This publication is to be used solely for the purpose of reference or operation of our software system. No part of this work can be disseminated for other purposes.

In the event of inadvertent or deliberate publication, we intend to enforce its right to this work

under copyright laws as a published work. Those having access to this work may not copy, use, or disclose the information in this work unless expressly authorized by our company.

All information contained in this publication is believed to be correct. We shall not be liable for

errors contained herein nor for incidental or consequential damages in connection with the furnishing, performance, or use of this material. The information this publication refers to is protected by copyrights or patents and does not convey any license under the patent rights of our company, nor the rights of others. We do not assume any liability arising out of any infringements of patents or other rights of third parties.

Content of the manual is subject to change without prior notice.

Model: MD300K2

Version: Ver 4.0K2

Revised Date: July 24, 2019

# CONTENT

<b>1 INTRODUCTION</b>	4
1.1 Brief Introduction	4
1.2 Intended Use	4
1.3 Measurement Principle	4
1.4 Safety Information	4
1.5 Electromagnetism Interference	6
1.6 Explanation of Symbols	7
1.7 Product Features	7
1.8 Contraindication	7
<b>2 GENERAL DESCRIPTION</b>	8
2.1 Appearance	8
2.2 Power Supply	9
<b>3 TIME AND DATE SETTINGS</b>	10
<b>4 TAKE A MEASUREMENT</b>	10
4.1 Install the Probe	10
4.2 Measurement Screens	11
<b>5 SETTINGS</b>	13
5.1 System Setting	13
5.2 Alarm Setting	13
<b>6 DATA MANAGEMENT</b>	15
6.1 Data Review	15
6.2 SpO <sub>2</sub> Trend	15
6.3 PR Trend	16
6.4 Delete all data	16
<b>7 ALARM</b>	17
<b>8 MAINTENANCE AND REPAIR</b>	18
8.1 Maintenance	18
8.2 Safety Checks	18
8.3 Troubleshooting	20
8.4 Warranty and Repair	20
<b>APPENDIX A SPECIFICATIONS</b>	23
<b>APPENDIX B CLINICAL STUDY SUMMARY</b>	26
<b>APPENDIX B DECALARATION</b>	28

## 1 INTRODUCTION

### 1.1 Brief Introduction

Thank you for purchasing the MD300K2 pulse oximeter. The main functions of the device include SpO<sub>2</sub>, PR and PI (Perfusion Index) measurements, visual and audible indication, probe off and finger out indication, data storage, transmitting and review etc. Please read this manual carefully before using the device.

#### **Notes:**

- The illustrations applied in the manual may differ slightly from the actual device.
- The device is designed of handheld structure, please be sure not to turn upside down when using it.

### 1.2 Intended Use


The MD300K2 Pulse Oximeter is intended for continuous monitoring, spot-checking of oxygen saturation (SpO<sub>2</sub>), pulse rate (PR) and PI (Perfusion Index) of single adult, adolescent, child and infant patients in hospitals and clinics.

### 1.3 Measurement Principle

Principle of the oximeter is as follows: A mathematical formula is established making use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin(RHb) and Oxyhemoglobin (HbO<sub>2</sub>) in red and near-infrared zones. Operation principle of the instrument: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning and Recording Technology, so that two beams of different wavelength of lights (660nm red and 905nm near infrared light) can be focused onto a human nail tip through a clamping finger-type sensor. A measured signal obtained by a photosensitive element, will be shown on the oximeter's display through process in electronic circuits and microprocessor.

#### **Diagram of Operation Principle**

1. Red and Infrared-ray Emission Tube
2. Red and Infrared-ray Receipt Tube

### 1.4 Safety Information

#### **Conception of Warning, Caution and Note**

The Warning, caution and Note at this document are special information in favor of users operation.

**Warning** - Information concerning something that could possibly hurt the patient or operator.

**Caution** - Reminds the user to pay close attention to device operation, failure of which may cause abnormal function of the instrument.

**Note** - Informs the user of other important information by suggestion, requirement and supplement.

#### **Warnings!**

1. Before use, carefully read the manual. This device is intended for use by persons trained in professional health care. Our company will assume no warranty for using this equipment improperly.
2. Operation of the handheld pulse oximeter may be affected by the use of an electrosurgical unit (ESU).

3. Sensor malfunction may cause inaccurate data possibly resulting in patient injury or death, so pay close attention to the sensor and inspect it often.
4. Do not use the handheld pulse oximeter in an MRI or CT environment.
5. Although the pulse oximeter has alarms, it is not suggested for long time continuous monitoring.
6. Do not use the handheld pulse oximeter in an explosive atmosphere.
7. The handheld pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
8. Check the pulse oximeter sensor application site every half an hour to determine the positioning of the sensor and circulation and skin sensitivity of the patient.
9. When begin to gauge associated body features please follow the doctor's counseling.
10. When link this equipment to other peripherals, make sure you are sophisticated operator to handle this device. Any peripherals should be in the light of protocol of IEC 60950 or IEC 60601-1. Any input/output device should be following the protocol of IEC 60601-1.
11. Do not sterilize the device using autoclaving, ethylene oxide sterilizing, or immersing the device in liquid. The device is not intended for sterilization.
12. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.
13. This equipment complies with IEC 60601-1-2:2014 for electromagnetic compatibility for medical electrical equipment and/or systems. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of such interference due to close proximity or strength of a source might disrupt the performance of this device.
14. You should operate the equipment according to the EMC information provided in the accompanying documents.
15. Portable and mobile RF communications equipment can affect medical electrical equipment.
16. This equipment should not be used adjacent to or stacked with other equipment.
17. This equipment is not intended for use during patient transport outside the healthcare facility
18. When connecting this device to other peripherals, make sure that you are qualified to operate this device. Any peripheral must be certified according to the protocol of IEC 60601-1. Any input/output device should follow the protocol of IEC 60601-1.
19. When using the equipment, the voice of the environment is not greater than 45 db.

**Rx only: "Caution: Federal law restricts this device to sale by or on the order of a physician."**

#### **Cautions:**

1. The pulse oximeter must be able to measure the pulse properly to obtain an accurate SpO<sub>2</sub> measurement. Verify that nothing is hindering the pulse measurement before relying on the SpO<sub>2</sub> measurement.
2. Worn-out data cables may also cause inaccurate data, so if the data is used as a reference to treat a patient, pay special attention to data cable and check it more frequently.
3. Do not tangle the SpO<sub>2</sub> cable with the wires of ES (Electrosurgery) equipment.

4. Single use accessories should never be re-used.
5. Only use SpO<sub>2</sub> sensors specified by the manufacturer. Other SpO<sub>2</sub> sensors may cause improper performance.
6. Unplug the sensor from the monitor before cleaning or disinfecting to prevent sensor or monitor from being damaged, and to prevent user under safety situation.
7. Alarm must be set up according to different situation of individual patient. Make sure that audio sound can be activated when alarm occurs.

**Notes:**

1. Optical cross talk can occur when two or more sensors are located in adjoining areas. It can be eliminated by covering each site with opaque material. Optical cross talk may adversely affect the accuracy of the SpO<sub>2</sub> readings.
2. Obstructions or dirt on the sensor's red light or detector may cause a sensor failure. Make sure there are no obstructions and the sensor is clean.
3. For routine equipment maintenance, please refer to the service procedures at the associated section as indicated in the manual.
4. As to the other concerns for attention, please carefully look through the specific chapter in this instruction.

**Inaccurate measurements may be caused by:**


1. Significant levels of dysfunctional hemoglobin (such as carbonyl - hemoglobin or methemoglobin);
2. Intravascular dyes such as indocyanine green or methylene blue;
3. High ambient light. Shield the sensor area if necessary;
4. Excessive patient movement;
5. High-frequency electrosurgical interference and defibrillators;
6. Venous pulsations;
7. Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line;
8. The patient has hypotension, severe vasoconstriction, severe anemia, or hypothermia;
9. The patient is in cardiac arrest or is in shock;
10. Fingernail polish or artificial fingernails;
11. Weak pulse quality (low perfusion);
12. Low hemoglobin;

**1.5 Electromagnetism Interference**

This oximeter is designed and tested in compliance with the EMC standard, complying with the international standard for the EMC of the electronic medical device – IEC 60601-1-2. However, because of the proliferation of radio frequency transmitting equipment and other sources of electrical noise in the health-care and home environments (e.g. cellular phones, mobile two-way radios, electrical appliances) it is possible that high levels of such interference due to close proximity or strength of a source, may result in disruption of performance of this device.

This apparatus complies with the IEC 60601-1-2 international standard. The requirements of this international standard are: CISPR11, GROU1, and CLASS B.

## 1.6 Explanation of Symbols

<i>Symbol</i>	<i>Explanation</i>	<i>Symbol</i>	<i>Explanation</i>
	Type BF applied part		Protected against dripping water
	Caution		Prevent from rain
	Date of Manufacture		Manufacturer's information
	Serial number		User ID
	Audio alarm inhibition		Battery power indication
	Storage temperature and relative humidity		Follow instructions for use
	Pulse rate		Hemoglobin Oxygen Saturation
	USB cable is connected		Battery cover unlock / lock
	European union approval		Authorized representative in the European community
	Do not discard the device and other components		Class II equipment
	Waste electrical and electronic equipment		Beep silence
	Indicate the signal is not stable		

## 1.7 Product Features

- Compact, light in weight and simple to use.
- TFT displaying screen with adjustable backlight.
- Up to 127 user ID and 72-hour records storage.
- Visual and three-level audible alarms, low battery alarm.
- Data transfer to PC by USB cable.
- Powered by three AA alkaline batteries or power adapter (optional).

## 1.8 Contraindication

None

## 2 GENERAL DESCRIPTION

The handheld pulse oximeter adopts 2.4 inch TFT displaying screen. It can display the SpO<sub>2</sub>, PR value, pulse bar, SpO<sub>2</sub> waveform and other kinds of parameters, such as the time, ID number, battery status and so on.

### 2.1 Appearance


Fig.2-1

#### **Description of Fig.2-1:**

1. Connector for probe and USB cable: To connect with SpO<sub>2</sub> probe and connect USB cable with PC for data transmission.
2. Alarm indicator: it will flash with yellow when alarms occur.
3. Displaying screen.
4. Navigation buttons (UP, DOWN, LEFT, RIGHT): the function of the left button is confirming the items that you choose; the function of the right button is back to the previous or menu interface. In the measuring interface, press the right button to shift between digital mode and wave mode.
5. Menu button: short press to enter the main menu; long press it to silence the audio alarm, the device will be appear the audio alarm indicator icon covered with "X" and countdown, long press it again, the device will recover the audio alarm.
6. Power button: Press and hold the button for 3 seconds to power the device on, and for about 4 seconds to power it off.
7. Connector for adapter.
8. USB interface: Used for update the software, only serves engineers;
9. Side UP/DOWN button: Short press it to increase/decrease the volume of alarm; Press it for extended time to increase/decrease the brightness of the backlight.

**NOTE:** In the measuring interface, press the menu button for extended time, the alarm pause will be turn on; press it again, the alarm pause will be turn off.

## 2.2 Power Supply

The device can be powered by 3 AA alkaline batteries.


Fig.2-2

### Batteries Installation:

- 1) **Open the battery cover:** Slide the fixing screw slightly in the rear panel to the position which is marked with “” and then open the cover.
- 2) **Batteries installation:** Install 3 AA alkaline batteries or rechargeable batteries lightly as indicated by the polarity signs as shown in Fig.2-2.
- 3) **Close the battery cover:** Close the battery cover and slide the screw to the  position. It means that the battery cover is locked.

### Battery life and replacement

When the battery voltage is lower, the battery indicator will become empty and its frame shows in red, which means few of battery capacity remains. You should replace the batteries with new ones in time.

### Warnings!

- If battery fluid should get into your eyes, immediately rinse with plenty of clean water. Consult a doctor immediately.
- Be sure to insert the batteries in the correct polarities, as indicated by polarity markings inside the battery compartment.

### Cautions!

- Do not use batteries not specified for this unit.
- Do not dispose of batteries in fire.
- Remove the batteries from this unit when you are not going to use it for a long period.
- Do not use batteries of different types together.
- Do not use new and used batteries together.
- Do not use the broken battery.
- Dispose of batteries in accordance with the local ordinances and regulations.

## AC Power Supply

The device can be supplied by AC power through connecting the device to AC adapter. Note: Use the AC power supply, make sure put the device in the safety and proper place and convenient to power off.

### Warnings!

- Be sure to use the adapter that specified for this device.
- Plug and unplug the adapter cautiously to avoid injuries caused to your body.
- If the device suddenly power off, please take out your finger at once, and then connect power or install the batteries.

## 3 TIME AND DATE SETTINGS

Set the correct time according to the following steps:

- 1) Press the power switch for 3 seconds to power on the oximeter and then press the menu button to enter the main menu, refer to the Fig.3.1.


Fig.3.1


Fig.3.2

- 2) Press the Navigation button to select "Date and Time" item, and then press the OK button to enter the time setup screen, refer to Fig.3.2. Press the Navigation button adjust the value, and then press the OK button to confirm the value. The date is displayed as the order of Year-Month-Day and Time of Hour-Minute (time show: 24-hour).


## 4 TAKE A MEASUREMENT

### 4.1 Install the Probe

1. Before use, check the pulse oximeter for mechanical damage.


2. Insert the SpO<sub>2</sub> probe into the probe socket (refer to Fig.4.1).  
If the SpO<sub>2</sub> probe is disconnected from the unit, a prompt "Probe Off" will be shown in the status bar in the top left of the display screen.

**Note:** The probes that specified for this device conform to the ISO80601-2-61 standard. The connector is also applied to uploading data to MedView software for reviewing. For detailed information, refer to the instruction manual of MedView software.

**Fig.4.1**

After finishing the time setting, plug your finger into the sensor as shown in the following illustration.


**Fig.4.2** Placement of the sensor


Select the suitable sensor in terms of type and dimension.

## 4.2 Measurement Screens

There are two display modes as shown in the following pictures.

**Fig.4.3** (Digital Mode)

**Fig.4.4** (Wave Mode)

### Description of measurement screens:

1. **Measuring:** The pulse oximeter is in the status of measuring.
  - It shows "Finger off" when there is no finger inserted or no signal is detected.
  - It shows "Probe off" when the sensor/probe is not connected to the pulse oximeter.
2. **SpO<sub>2</sub>%:** SpO<sub>2</sub> area of display
  - It shows the oxygen saturation level of functional hemoglobin during normal measurement.
  - The color of the SpO<sub>2</sub> value will become red when the SpO<sub>2</sub> is beyond the alarm limits.
  - It shows two dashes throughout probe off and finger out conditions.
3. **100:** SpO<sub>2</sub> high alarm limit; **90:** SpO<sub>2</sub> low alarm limit.

4. :  Pulse bar
5. **100**: PR high alarm limit; **060**: PR low alarm limit.
6. **PR**: PR area of display
  - It shows the pulse rate in beats per minute during normal measurement.
  - The color of the PR value will become red when the PR is beyond the alarm limits.
  - It shows three dashes throughout probe off and finger out conditions.
7. **ID**: 001, the ID number of the current patient is 001.
8. 12:00: The current time.
9. **PI%**: Perfusion Index indicator display area.

### Warnings!

- Use only SpO<sub>2</sub> probes provided by manufacturer. Other SpO<sub>2</sub> probes may cause improper performance.
- Do not use an SpO<sub>2</sub> probe with exposed optical components.
- Excessive patient movement may cause inaccurate measurements.
- Tissue damage can be caused by incorrect operation or misusing probe; for example, by wrapping the probe too tight. Inspect the probe site to ensure the skin's integrity and the adhesion position of the probe is correct. More frequent inspection should be taken if necessary.
- Loss of pulse signal can occur in any of the following situations:
  - a) The probe is too tight;
  - b) There is excessive illumination from light sources such as a surgical lamp, a bilirubin lamp, or sunlight; Pulse sensor should obviate the light source, e.g. radial lamp or infrared lamp.
  - c) A blood pressure cuff is inflated on the same extremity as the one to which an SpO<sub>2</sub> sensor is attached.
- After measurement, please take off your finger and press and hold the power button to turn off the device.

## 5 SETTINGS

### 5.1 System Setting

Pick and enter the [System Setting] interface from the main menu. And then press the Navigation buttons to select different item to set.


Fig.5.1


Fig.5.2

**Alarm Volume:** you can adjust the value of alarm volume, there are 7 levels, and the default level is 3

**Alarm Pause:** there are two modes, 60s and 120s, and the default mode is 120s. The device audibly alarm upon new alarm conditions.

**Beep Tone:** the level is from 0 to 7, and the default level is 3.

**Backlight Setting:** the level of brightness is from 1 to 7, and the default level is 3.

**Language:** English, French, German, Spanish, Italian, Japanese, Russian and Chinese.

**Screen Sleep Mode:** 1minute, 10 minutes, 30 minutes, screen always on, and the default is 1 minutes

**Charging management:** Charging activated, Charging stop.

**ID setup:** ① press OK button, ② press Navigation button to change the number, ③ press OK button to confirm.

**Factory Default:** recover to factory reset.

#### Note:

1. Every time enter into the Alarm Volume and Alarm Pause, you should input the password, refer to Alarm Setting.
2. The ID range is 1~127.

### 5.2 Alarm Setting

Pick the [Alarm Setting] interface from the main menu. Before setting, please enter password (1234) to set the parameter. Or you can direct access to check the parameter but not to change.

## How do you input passwords?

1. Press the Navigation button to change the numbers.
2. Press the OK button to confirm the number.
3. Press the Navigation button switch to the next numbers.
4. Press the OK button to confirm your selection.
5. Repeat the step one.

## How to change the passwords?

In the measuring interface, hold and press the menu button for 5s to change the passwords.

The first time, enter the old passwords.

The second time, enter the new passwords.

Then you can change the passwords.

### Notes:

1. The read-only password is 0000. Under this password, you can only check the parameter but not change.
2. The make-changes password is 1234, enter this password, you can set the parameters.
3. You can change the password. If you forget the password, Please choose "Factory Default" in System Setting, the password will Recover to factory password (1234) .
4. Every time enter into the Alarm Setting in the main menu, Alarm Volume and Alarm Pause in the System Setting, you should input the password.

Enter the [Alarm Setting] interface.


Fig.5.4

You can select the Alarm on or off.

High limit SpO<sub>2</sub> range is 71~100, Low Limit SpO<sub>2</sub> range is 70~99

High Limit PR range is 31~250, Low Limit PR range is 30~249


Fig.5.3

## 6 DATA MANAGEMENT

Pick and enter the [Data Management] interface from the main menu. You can choose continuous measurement or not.


Fig.6.1

### 6.1 Data Review

Pick and enter the “Data review” interface as shown in Fig.6.1. By pressing the UP / DOWN button, the user can turn to the earlier/later item page.

The pulse oximeter can record the alarming parameter marked with red color. Press the MENU button, it will return to the main menu interface.

The screenshot shows a device interface with a status bar at the top displaying 'Measuring', a signal strength icon, '120s', and a battery level icon. Below the status bar, there is a table with the following data:

Time	SpO <sub>2</sub>	PR	ID
23/04 06:00:20	98	70	1
23/04 06:00:16	98	70	1
23/04 06:00:12	98	70	1
23/04 06:00:08	98	70	1
23/04 06:00:04	90	60	1
23/04 06:00:00	90	60	1
23/04 05:59:56	90	60	1
23/04 05:59:52	90	60	1
23/04 05:59:48	90	60	1
23/04 05:59:44	90	60	1

At the bottom of the table, there are navigation buttons: 'Page 01/80' and 'Back'.

Fig.6.2

### 6.2 SpO<sub>2</sub> Trend

Pick and enter the “SpO<sub>2</sub> Trend” interface as shown in Fig.6.3. Press the Left / Right button; you can turn to the earlier or later item page. Press the Up button, the pulse oximeter returns to the previous interface. Press the MENU button, the pulse oximeter returns to the main menu interface.


Fig.6.3

On the above of the trend, the date and time of the first item are displayed, with month/day; hour: minute: second.

### 6.3 PR Trend

Pick and enter the “PR Trend” interface as shown in Fig.6.4. Press the Left / Right button; you can turn to the earlier or later item page. Press the Up button, the pulse oximeter returns to the previous interface. Press the MENU button, the pulse oximeter returns to the main menu interface.


Fig.6.4


### 6.4 Delete all data

Pick and enter the “Delete all data” interface as shown in Fig.6.5. You can select “Yes” or “No” by pressing the up or down button, and by pressing the OK button to confirm your selection.

**Note:**

1. Please take caution to the deletion of data; you will never get the data back once deleted.
2. Data is unable to be deleted during measurement.

**Fig.6.5**

## 7 ALARM

### ALARM PRIORITY:

There are two-level priorities for selection.

High priority: the highest level alarm, indicates the patient is in the very dangerous situation.

Low priority: indicate the technical alarm caused by the device itself.

Alarms of the oximeter include technical and physiological alarms. All the two priorities divided by built-in module and cannot be changed by user.

### Assignment of priority:

	High	Low
<b>Paramter</b>	SpO <sub>2</sub> /PR	/
<b>Value</b>	Red	/
<b>Alarm lamp</b>	Flashing	/
<b>Lamp Frequency</b>	1.5Hz	/
<b>Audiblesound</b>	Di- Di – Di ----- Di - Di	Di
<b>Alarm cycle</b>	3 s	20 s
<b>Alarm info</b>	SpO <sub>2</sub> too high/low, PR too high/low, Battery power low	Probe off/Finger off

### Notes:

1. The alarm will appear if the measurement value out of range.
2. The alarm sound will go on until alarm disappears or is turned off.
3. After silencing the alarm, the corresponding indicator will indicate this.
4. The power low alarm: the corresponding indication lamp will be flashing with a red frame.

## Warnings!

***When the alarm occurs, check patients' conditions immediately.***

- When the alarm occurs, check patients' conditions immediately.
- Check which parameter is alarming or which alarm is going on.
- Check patient's condition.
- Search for the source of alarm.
- Make the alarm mute if necessary.
- Check the alarm when no warning.

### **Alarm delay**

The alarm condition delay and alarm signal generation delay: less than 1s.

#### **Note:**

1. The pulse rate correspondence with the user's pulse rate. It based on the user's actual pulse rate.
2. Use the alarm setting in different areas will be brings the potential dangerous.
3. The alarm setting can recover if the power-break time less than 30s.
4. Set the high parameter value with simulator to test the efficiency of the alarm system.
5. Do not set the parameter value out of the range, or the alarm system will failure.
6. The device can reserve alarm setting if it power break.

## **8 MAINTENANCE AND REPAIR**

### Warnings!

- The advanced circuit inside the oximeter does not require periodic calibration and maintenance, except replacing the batteries.
- Don't open the cover of oximeter or repair electronic circuits. Its open will cause the damage of the device and the annulment of the guarantee.

### **8.1 Maintenance**

Use only the substances approved by us and methods listed in this chapter to clean or disinfect your equipment. Warranty does not cover damage caused by unapproved substances or methods.

We make no claims regarding the efficacy of the listed chemicals or methods as a means for controlling infection. For the method to control infection, consult your hospital's Infection Control Officer or Epidemiologist. Keep your equipment and accessories free of dust and dirt. To avoid damage to the equipment, follow these rules:

- Always dilute according the manufacturer's instructions or use lowest possible concentration.
- Do not immerse part of the equipment into liquid.

- Do not pour liquid onto the equipment or accessories.
- Do not allow liquid to enter the case.
- Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).

**Caution:** If you spill liquid on the equipment or accessories, contact us or your service personnel.

**Note:** To clean or disinfect reusable accessories, refer to the instructions delivered with the accessories.

## 8.2 Safety Checks

Before every use, or after your pulse oximeter has been used for 6 to 12 months, or whenever your pulse oximeter is repaired or upgraded, a thorough inspection should be performed by qualified service personnel to ensure the reliability. Follow these guidelines when inspecting the equipment:

- Make sure that the environment and power supply meet the requirements.
- Inspect the equipment and its accessories for mechanical damage.
- Make sure that only specified accessories are applied.
- Inspect if the alarm system functions correctly.
- Make sure that the batteries meet the performance requirements.
- Make sure that the pulse oximeter is in good working condition.

In case of any damage or abnormality, do not use the pulse oximeter. Contact your hospital's biomedical engineers or your service personnel immediately.

### Cleaning

Please use medical alcohol to clean the silicone touching the finger inside of SpO<sub>2</sub> probe with a soft cloth dampened with 70% isopropyl alcohol. Also clean the being tested finger using alcohol before and after each test. To clean your equipment, follow these rules:

1. Shut down the pulse oximeter and take the batteries out of the battery wharf.
2. Clean the display screen using a soft, clean cloth dampened with a glass cleaner.
3. Clean the exterior surface of the equipment using a soft cloth dampened with the cleaner.
4. Wipe off all the cleaning solution with a dry cloth after cleaning if necessary.
5. Dry your equipment in a ventilated, cool place.

### Disinfecting

The applied parts touching the patients' body are required to be disinfected once after each use. The recommended disinfectants include: ethanol 70%, isopropanol 70%, glutaraldehyde 2% liquid disinfectants.

Disinfection may cause damage to the equipment and is therefore not recommended for this pulse oximeter unless otherwise indicated in your hospital's servicing schedule. Clean the pulse oximeter before disinfecting it.

**Caution:** Never use EtO or formaldehyde for disinfection.

"dudu" voice should be heard.

## 8.3 Troubleshooting

<b>Malfunction</b>	<b>Reason</b>	<b>Solution</b>
SpO <sub>2</sub> or PR is shown unstably	<ol style="list-style-type: none"><li>1. Finger might not be plugged deep enough.</li><li>2. Finger is trembling or patient's body is in movement status.</li></ol>	<ol style="list-style-type: none"><li>1. Retry by plugging the finger.</li><li>2. Try not to move.</li></ol>
The Oximeter cannot be powered on	<ol style="list-style-type: none"><li>1. Power of batteries might be inadequate or not be there at all.</li><li>2. Batteries might be installed incorrectly.</li><li>3. The Oximeter might be damaged.</li></ol>	<ol style="list-style-type: none"><li>1. Please replace batteries</li><li>2. Please reinstall the batteries</li><li>3. Please contact with local customer service centre.</li></ol>
Indication lamps are suddenly off	<ol style="list-style-type: none"><li>1. The product is automatically powered off when no signal is detected longer than 8 seconds</li><li>2. Lower power</li></ol>	<ol style="list-style-type: none"><li>1. Normal</li><li>2. Replace the batteries</li></ol>
Others		Please contact with local customer service centre.

## 8.4 Warranty and Repair

### 8.4.1 Maintenance Method

Repair response time: AM9:00 to PM17:30 on Monday to Friday except legal holiday.

Repair time: AM9:00 to PM17:30 on Monday to Friday except legal holiday.

### 8.4.2 Exempt and limitation

a) Our company isn't responsible for such damage caused by force nature. For example: fire, thunder flash, flood, cyclone, hail, earthquake, house collapse, commotion, plane failing and traffic accident, deliberate damage, lack of fuel or water, labor and capital bother, strike and stop-working etc.

b) No-service offer

- The cost and insurance charge of disassembling, refurbishing, repackaging and conveying of the oximeter or the part of it.
- Damage or loss sustained due to inspected or repaired by other institute that is not certified.
- The damage and failure caused by user or its representative who doesn't use the device according to the operator's manual

c) The damage or loss sustained due to connection to peripheral equipment (such as printer, computer etc.), that are not provided by our company are not covered by the warranty.

d) Responsibility limitation

In the duration of warranty, if user changes the parts manufactured by other manufacturers without our company permission, our company is entitled to cancel contract.

### **8.4.3 User Guarantees**

a) Please read user manual carefully before operation

b) Please operate and make daily maintenance as request of manual and guarantee

c) Power supply and environment must be maintained under manual specifications.

### **8.4.4 No-guarantee principle**

- The device does not remain in original condition.
- The shell of the device is breached or cracked.
- Evidence of water damage.
- Accessories adulterated or appearance of physical abuse.
- Evidence of crushing damage to the probe.
- Original Packaging during transportation is not used.
- Non authorized service is performed on oximeter.
- Damage to a product as a result of not conforming to manual specifications.
- The work environment is not eligible.
- There is smear or marks that are not belong to the instrument and cannot be removed from the outside surface of the instrument.
- The circuit is short and damaged due to liquor or other stuff flow in the instrument or its fittings.
- All probe and its accessories are not free replacement.
- If any code label of parts is damaged or missing, this warranty shall become null and void.
- For example of code label.

Such damage of probe caused by mechanical force doesn't belong to free change range.

- During measurement of SpO<sub>2</sub>, principle leads to measuring value difficultly or inaccurate measurement.
- Maintenance seal of oximeter are not opened.

### **8.4.5 User's Special Request for Guarantee Time**

As we stipulate the warranty period according to the relevant electronic regulation of country, which we stipulate is one year, accessory is three months. When customer requires to extend the warranty period, you should consider whether it is reasonable. Because electronic product is quickly replaced, as to the warranty period over three years, purchased accessories may be out of stock. In this case, we will adopt to entirely upgrade or replace the old, you should pay the minimum acceptable cost of renewed device.

### **8.4.6 Repackage**

- Take all the accessories and put them into plastic cover

- Try to use original package and packing material. User will be responsible for such damage caused by bad package during transportation.
- Please offer guarantee list and copy of invoice to standby with the period of guarantee.
- Please describe failure phenomenon in detail and altogether offer oximeter.

### ***Storage and Transportation***

Storage: Storage Temperature  $-20^{\circ}\text{C}\sim 70^{\circ}\text{C}$ , Relative Humidity  $\leq 93\%$

Transportation: Transport by airline, train or vessel after packing according to request.

Package: We pack the product with the hard bag. We put the foam between the inner box and the cartoon to alleviate the shake.

## APPENDIX A SPECIFICATIONS

### Notes:

- Specifications may be changed without prior notice.
- The circuit diagrams, the list of components, the illustration of diagrams, and the detailed rules of calibration, are provided exclusively to professional personnel authorized by our company.
- The equipment has been calibrated, users do not need to calibrate. In order to ensure the accuracy of the probe, please change the probe once a year. Make sure that the type of probe need to be specified.

### Display

Data: SpO<sub>2</sub>, PR, pulse bar, PI, pulse waveform

Others: Information such as connection status of probe and power low indication.

### Alarm

Alarm: SpO<sub>2</sub> and pulse rate value, probe off, battery exhausted

Alarm mode: audio alarm, visual alarm and information

Alarm limits range: SpO<sub>2</sub> 70%~100%, PR 30bpm~250bpm

Default limits: SpO<sub>2</sub> High 100%, low 95%; PR High 100bpm; low 60bpm

### SpO<sub>2</sub>

Measurement range: 70%~100%

Resolution: 1%

Accuracy:

Classification		Accuracy
No motion	Adults	±2%
	Pediatrics	±2%
	Neonates	±3%
Motion	Adults	±3%
	Pediatrics	±3%
	Neonates	±3%
Low Perfusion	Adults	±2%
	Pediatrics	±2%
	Neonates	±3%
Low Perfusion Performance	PI=0.1%	±3%
Anti-interference Performance	Fluke Index2 simulator Level2	±3%

## Pulse Rate

Measurement range: 30~250 bpm

Resolution: 1 bpm

Accuracy:  $\pm 2$  bpm or 2% (The larger are applied)

## Perfusion Index

Display range: 0.1%~20%


Accuracy: 0.1%~1.0%,  $\pm 0.2\%$ ; 1.1%~20%,  $\pm 20\%$

## Probe LED Specifications

	Wavelength	Radiant Power
RED	660 $\pm$ 3nm	3.2mW
IR	905 $\pm$ 10nm	2.4mW

## Equipment data update period

As shown in the following figure. Data update period of slower average is 8s.

## Operation Environment

Operating temperature: 0°C ~ 40°C

Relative humidity:  $\leq 80\%$ , no condensation

Atmosphere pressure: 86 kPa ~ 106 kPa

## Transport and Storage Environment

Storage Temperature: -20°C ~ 70°C

Relative humidity:  $\leq 93\%$ , no condensation

Atmosphere pressure: 50 ~ 106 kPa

## Power Requirement

Power supply: Three AA alkaline or Ni-MH batteries, or power adapter

Working time: work for 10 hours continuously at most

## **AC Adapter**

Input Voltage: AC 100~240V

Input Frequency: 50Hz~60Hz

Output Voltage: DC 5V±5%

Output Current: 2A MAX

## **SMD chip fuse**

32V3A0466003.NR(1206)

## **Data store and replay**

Store and replay 72 hours SpO<sub>2</sub> and PR value, the time interval is 4 seconds.

## **Outline of product**

Dimension: 143.3mm(L)X67.4mm(W)X36mm(H)

Weight: ≤140±10g (excluding the batteries)

## **Box contents:**

1. Three AA alkaline batteries
2. One instruction manual
3. One adult finger probe: M-50E012CS09
4. One pediatric finger probe: M-50B008CS09 (optional)
5. One disposable probe: M-50J033CS045 (optional)
6. One USB cable
7. Adapter(optional)

## **Equipment Classification**

Classification according to IEC-60601-1:

**According to the type of protection against Electrical shock:** Internal electrical power source equipment and class II equipment;

**According to the degree of protection against Electrical shock:** Type BF equipment, applied part is SpO<sub>2</sub> probe;

**According to the degree of protection against harmful ingress of water:** IPX1;

**According to the methods of sterilization or disinfection:** Non-sterilization, use of liquid surface disinfectants only;

**According to the mode of operation:** Continuous operation.

## APPENDIX B CLINICAL STUDY SUMMARY

The following details are provided to disclose actual performance observed in the clinical validation study of healthy adult volunteers. The ARMS value analysis statement and Bland-Altman plot of data for MD300K2 and its supporting probes are shown as following:


### Clinical study details of MD300K2 Pulse Oximeter and its supporting M-50E012CS09

Oximeter probe:

#### ARMS Value Analysis Statement

Item	90--100	80--<90	70--<80
#pts	78	74	66
Bias	-0.73	-0.59	0.45
ARMS	1.46	1.80	1.99

#### Bland-Altman Plot Graphic


### Clinical study details of MD300K2 Pulse Oximeter and its supporting M-50B008CS09

Oximeter probe:

#### ARMS Value Analysis Statement

Item	90--100	80--<90	70--<80
#pts	78	74	66
Bias	-0.10	-0.31	-0.03
ARMS	1.19	1.40	1.82

## Bland-Altman Plot Graphic


Clinical study details of MD300K2 Pulse Oximeter and its supporting M-50J033CS045

Oximeter probe:

ARMS Value Analysis Statement

Item	90--100	80--<90	70--<80
#pts	78	74	66
Bias	-0.51	0.41	1.56
ARMS	1.34	1.49	2.36

## Bland-Altman Plot Graphic

## APPENDIX C DECLARATION

### Guidance and manufacturer's declaration - Electromagnetic emission---- for all EQUIPMENT AND SYSTEM

1	Guidance and manufacturer's declaration - electromagnetic emission		
2	The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic specified below. The customer or the user of the model MD300K2 Pulse Oximeter should assure that it is such an environment.		
3	Emissions test	Compliance	Electromagnetic environment-guidane
4	RF emissions CISPR11	Group 1	The model MD300K2 Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
5	RF emissions CISPR11	Group B	The MD300K2 Pulse Oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
6	Harmonic emissions IEC 61000-3-2	Class A	
7	Voltage fluctuations/ IEC 61000-3-3	Complies	

### Guidance and manufacturer's declaration - Electromagnetic Immunity- For all Equipment and Systems

Guidance and manufacturer's declaration- electromagnetic immunity

The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the model MD300K2 Pulse Oximeter should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-4	$\pm 2,4,6,8$ kV contact $\pm 2,4,8,15$ kV air	$\pm 2,4,6,8$ kV contact $\pm 2,4,8,15$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1$ kV differential mode $\pm 2$ kV common mode	$\pm 1$ kV differential mode $\pm 2$ kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-1	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 2 5 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 2 5 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model MD300K2 Pulse Oximeter requires continued operation during power main Interruptions, it is recommended that the model MD300K2 Pulse Oximeter be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

**Guidance and manufacturer's declaration- electromagnetic immunity-**  
**For EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING**

Guidance and manufacturer's declaration- electromagnetic immunity			
The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the model MD300K2 Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80kHz	3V	Portable and mobile RF communications equipment should be used no closer to any part of the model MD300K2 Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2\sqrt{P}$ 150kHz to 80 MHz $d=1.2\sqrt{P}$ 80MHz to 800MHz $d=2.3\sqrt{P}$ 800MHz to 2.5GHz
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	Where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacture and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE1 At 80MHz and 800MHz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

- a.) Field strengths from fixed transmitters, such as base situation for radio (cellular/cordless) telephones and land/mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model MD300K2 Pulse Oximeter is used exceeds the applicable RF compliance level above, the Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, the additional measures may be necessary, such as reorienting or relocating the model MD300K2 Pulse Oximeter.
- b.) Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM-for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING**

Recommended separation distances between portable and mobile RF communications equipment and the MD300K2 Pulse Oximeter

The model MD300K2 Pulse Oximeter is intended for use in an electromagnetic environment in which radiated RF disturbance is controlled. The customer or the user of the model MD300K2 Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the model MD300K2 Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter(W)	Separation distance according to frequency of transmitter (m)		
	150KHz to 80 MHz $d=1.2\sqrt{P}$	80MHz to 800 MHz $d=1.2\sqrt{P}$	800MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic interference is affected by absorption and reflection from structures, objects and people.

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-4	$\pm 2,4,6,8\text{kV}$ contact $\pm 2,4,8,15\text{kV}$ air	$\pm 2,4,6,8\text{kV}$ contact $\pm 2,4,8,15\text{kV}$ air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient/burst IEC 61000-4-4	$\pm 2\text{kV}$ for power supply lines $\pm 1\text{kV}$ for input/output lines	$\pm 2\text{kV}$ for power supply lines $\pm 1\text{kV}$ for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1\text{kV}$ differential mode $\pm 2\text{kV}$ common mode	$\pm 1\text{kV}$ differential mode $\pm 2\text{kV}$ common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-1	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) For 2.5 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) For 2.5 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model MD300K2 Pulse Oximeter requires continued operation during power main interruptions, it is recommended that the model MD300K2 Pulse Oximeter be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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## Guidance and manufacturer's declaration- electromagnetic immunity-

### For EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration- electromagnetic immunity			
The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the model MD300K2 Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80kHz	3V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model MD300K2 Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d=1.2\sqrt{P} \quad 150\text{kHz to } 80\text{MHz}$ $d=1.2\sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	<p><math>d=2.3\sqrt{P}</math> 800MHz to 2.5GHz</p> <p>Where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacture and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE1 At 80MHz and 800MHz, the higher frequency range applies.</p> <p>NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-4	±2,4,6,8kV contact ±2,4,8,15kV air	±2,4,6,8kV contact ±2,4,8,15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	±2kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-1	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 2.5 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 2.5 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model MD300K2 Pulse Oximeter requires continued operation during power main interruptions, it is recommended that the model MD300K2 Pulse Oximeter be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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**Guidance and manufacturer's declaration- electromagnetic immunity-  
For EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING**

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## 1.6 Symbole

<i>Symbol</i>	<i>Explanation</i>	<i>Symbol</i>	<i>Explanation</i>
	Type BF applied part		Protected against dripping water
	Caution		Prevent from rain
	Date of Manufacture		Manufacturer's information
	Serial number		User ID
	Audio alarm inhibition		Battery power indication
	Storage temperature and relative humidity		Follow instructions for use
<b>bpm</b>	Pulse rate	<b>SpO<sub>2</sub></b>	Hemoglobin Oxygen Saturation
	USB cable is connected		Battery cover unlock / lock
	European union approval		Authorized representative in the European community
	Do not discard the device and other components		Class II equipment
	Waste electrical and electronic equipment		Beep silence
<b>?</b>	Indicate the signal is not stable		

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-4	$\pm 2,4,6,8$ kV contact $\pm 2,4,8,15$ kV air	$\pm 2,4,6,8$ kV contact $\pm 2,4,8,15$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1$ kV differential mode $\pm 2$ kV common mode	$\pm 1$ kV differential mode $\pm 2$ kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-1	<5% UT (>95%dip in UT) for 0.5 cycles 40% UT (60%dip in UT) for 5 cycles 70% UT (30%dip in UT) For 2 5 cycles <5% UT (>95%dip in UT) for 5 sec	<5% UT (>95%dip in UT) for 0.5 cycles 40% UT (60%dip in UT) for 5 cycles 70% UT (30%dip in UT) For 2 5 cycles <5% UT (>95%dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model MD300K2 Pulse Oximeter requires continued operation during power main Interruptions, it is recommended that the model MD300K2 Pulse Oximeter be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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Guidance and manufacturer's declaration- electromagnetic immunity			
The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the model MD300K2 Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80kHz	3V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model MD300K2 Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d=1.2\sqrt{P} \quad 150\text{kHz to } 80\text{MHz}$ $d=1.2\sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	<p><math>d=2.3\sqrt{P}</math> 800MHz to 2.5GHz</p> <p>Where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacture and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE1 At 80MHz and 800MHz, the higher frequency range applies.</p> <p>NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-4	$\pm 2,4,6,8\text{kV}$ contact $\pm 2,4,8,15\text{kV}$ air	$\pm 2,4,6,8\text{kV}$ contact $\pm 2,4,8,15\text{kV}$ air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic transient/burst IEC 61000-4-4	$\pm 2\text{kV}$ for power supply lines $\pm 1\text{kV}$ for input/output lines	$\pm 2\text{kV}$ for power supply lines $\pm 1\text{kV}$ for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1\text{kV}$ differential mode $\pm 2\text{kV}$ common mode	$\pm 1\text{kV}$ differential mode $\pm 2\text{kV}$ common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-1	$<5\%$ UT ( $>95\%$ dip in UT) for 0.5 cycles $40\%$ UT ( $60\%$ dip in UT) for 5 cycles $70\%$ UT ( $30\%$ dip in UT) For 2 5 cycles $<5\%$ UT ( $>95\%$ dip in UT) for 5 sec	$<5\%$ UT ( $>95\%$ dip in UT) for 0.5 cycles $40\%$ UT ( $60\%$ dip in UT) for 5 cycles $70\%$ UT ( $30\%$ dip in UT) For 2 5 cycles $<5\%$ UT ( $>95\%$ dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model MD300K2 Pulse Oximeter requires continued operation during power main Interruptions, it is recommended that the model MD300K2 Pulse Oximeter be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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Guidance and manufacturer's declaration- electromagnetic immunity			
The model MD300K2 Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the model MD300K2 Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80kHz	3V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the model MD300K2 Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d=1.2\sqrt{P}$ <p>150kHz to 80 MHz</p> $d=1.2\sqrt{P}$ <p>80MHz to 800MHz</p> $d=2.3\sqrt{P}$ <p>800MHz to 2.5GHz</p>
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	<p>Where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacture and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE1 At 80MHz and 800MHz, the higher frequency range applies.            NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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