

Product Data Sheet

Powderfree Nitrile, 3mil., finger textured, Green Tea

Material:	100% Nitrile (Acetonitrile-butadiene)	Revision:	1
Appearance:	Powder-free, non-sterile	Date:	06.04.2020
Thickness & Color:	4.0g Veggie Green (3.0mil)		
Scent:	Green Tea		
Surface:	Finger textured		

Dimension

Size	Palm Width (mm)	Lenth (mm)	Item Code
Extra-Small	70-80	240 (min.)	9883760
Small	81-90	240 (min.)	9883761
Medium	91-100	240 (min.)	9883762
Large	101-110	240 (min.)	9883763
Extra-Large	111-120	240 (min.)	9883764

Thickness

Point of Thickness Measurements	Single Wall (mm)
Finger (at 15 mm from the extreme tip)	0.10 – 0.16 mm
Palm (at the center of palm)	0.07 – 0.11 mm

Physical Properties

Criteria	Before Aging	After Aging
Tensile Strength (MPa)	≥ 14	≥ 14
Elongation (%)	500 (min.)	400 (min.)
Force at Break (N)	≥ 6	≥ 6

Powder Residue: Contains less than 2 mg per glove.
 Design Feature: Ambidextrous, straight fingers, finger textured and beaded cuff.
 Product Conformance: Conform to ASTM D6319, EN 455 and AS/NZS 4011 I
 In compliance with Medical Device Regulation MDR 2017/745 (CE class I)
 Packaging: 100 gloves per dispenser box and 10 dispensers per carton.
 Glove shelf-life: 3 year from manufacturing date

Pre-shipment

Inspection Criteria	Inspection Level	AQL
Dimension	13	Median Value must Comply-
Physical Properties	13	Median Value must Comply-
1000 mL Water leak Test	G-1	1.5
Visual Inspection or Major Defects	G1	2.5
Visual Inspection or Minor Defects	G1	4.0

Weight (for reference only)

Size	Weight
Medium	4.0 ± 0.2g

Quality Assurance: Manufacturing process in compliance with ISO 13485 Quality

Note: The data in above labels represent typical analyses of the glove at time of publication and not to be used as actual product specification.